

# Collection, Concentration and Long Term Room Temperature Stabilization of Forensic DNA in Liquid Format is a Reality

Geoffrey Routh, PhD.<sup>1</sup> & Bassam El-Fahmawi, PhD.<sup>2</sup>

1. BioChain Institute Inc. Newark California, USA, 94560. 2. Mawi DNA Technologies, Hayward, California, USA, 94545



## 1. Introduction

Crime scene investigators are trained to air-dry evidentiary samples following swab-based collection. Drying swabs can lead to irreversible binding of DNA to the swab material resulting in low DNA recovery. However this step can be challenging, especially in areas of high humidity and air moisture content. There are several factors affecting DNA integrity and yields of evidentiary samples:


- Irreversible binding to swab material
- Nucleases
- Microbial contamination
- Temperature fluctuations
- Heat
- Humidity
- UV light
- Improper collection
- Storage conditions and length of backlog
- Transport conditions & transit time

Mawi has developed iSWAB-ID, an efficient liquid-based sample collection system which utilizes swabs. iSWAB-ID enables long term room temperature stabilization of the collected sample at the point of collection while ensuring proper chain of custody. This system maximizes sample recovery and simplifies sample processing extensively over current practices, allowing for Enhanced ID Profiling.

## 4. Results

### I- DNA Extraction and Human DNA Quantification iSWAB-ID Efficiently Recovers and Stabilizes Human DNA from Both Reference & Evidentiary Samples in Long Term, Room Temperature Storage from Low, Medium and High Copy Number Samples

DNA concentration from swabs post iSWAB-ID collection: ND



## Evidentiary Material Processing Bottlenecks

Step #	Description	Time Spent (mins)
1	Screen	
2	Identify	
3	Collect	
4	<b>Dry</b>	<b>30 min/Sample</b>
5	Package	
6	Reporting and documentation	
7	Transport	
8	Store	
9	<b>Cut swab/ Punch FTA Card</b>	<b>10 min /Sample</b>
10	Lyse	
11	Extract	
12	Quantify	
13	<b>STR profiling</b>	<b>20-60% Further Analysis</b>
14	Data Analysis	
15	Report	

## 2. Objectives

- To assess the efficiency of iSWAB-ID in the collection, concentration, and long-term room temperature stabilization of both reference and evidentiary samples by assessing the usability of purified DNA in human ID-STR profiling assays.
- To assess the stability of iSWAB-ID samples during multistop-global transport process

## 3. Experimental Design

Collect reference and mocked evidentiary samples into the iSWAB-ID Sample Collection Device (400 µL)

↓


Transport to processing lab by standard mail

↓


Store samples at room temperature

↓


100 µL aliquoted from each iSWAB-ID sample post collection



### Lipstick on a Cup Rim: An Example of Touch DNA Sample



### Simulated Rape



## 5. Summary & Conclusions

### Evidentiary & Reference Samples Processing Bottlenecks: All Resolved with iSWAB-ID

Step #	Description	Time Spent (mins)
1	Screen	
2	Identify	
3	Collect	
4	<b>No Drying Required</b>	<b>0 min/Sample</b>
5	Package	
6	Reporting and documentation	
7	Transport	
8	Store	
9	<b>Cut swab/ Punch FTA Card</b>	<b>0 min/Sample Not Applicable</b>
10	Lyse	
11	Extract	
12	Quantify	
13	<b>STR profiling</b>	<b>5-10% Further Analysis</b>
14	Data Analysis	
15	Report	

Significant Improvement on Processing Efficiency:  
 ✓ Faster Collection  
 ✓ Higher First Pass Rate  
 ✓ Higher Sample processing Throughput

- The iSWAB-ID sample collection device efficiently recovered and stabilized DNA of forensic significance in liquid format at the point of collection.
- DNA collected and stabilized in iSWAB-ID at ambient temperature remained of sufficient quality to analyze for at least 6 months.
- Unlike samples processed from swabs, samples collected with iSWAB-ID yield more DNA allowing multiple analytical runs and sufficient material for archiving.
- Extended ambient temperature shelf-life allows cost and space-saving storage and eliminates sample degradation resulting from excessive backlogs.